

Rosé

Geriatric Chair



Rosé Geriatric Chair



A comfortable and versatile nursing aid

- **comfortable**
- **versatile**
- **easy-to-manoeuvre**
- **wide selection of accessories**

Rosé Geriatric Chair is comfortable and functional. It has been especially designed for use in nursing of the elderly and for patients requiring long-term care. It features a contoured seat which provides sturdy support and optimum seating comfort during extended use. Dual, swivelling castors make the chair easy to manoeuvre. The Rosé Geriatric Chair comes with several accessories, all of which are designed to make nursing tasks easier.

TECHNICAL SPECIFICATIONS

- Frame of epoxy-coated, formed steel tube
- Detachable padding, easy to keep clean, colour blue or beige
- Adjustable seating position
- Adjustable seating angle for easy access
- Dual, swivelling castors Ø 100 mm, rear castors equipped with brakes and steering facilities
- Separately liftable armrests
- Adjustable, contoured neck cushion
- Turn-up foot rest

DIMENSIONS

- Max. width 650 mm
- Seat width 490 mm
- Max. height 1250 mm
- Seat height 420 mm
- Max. length 750 mm
- Seat length 520 mm

ACCESSORIES

- separate dining table (pictures 3 and 4)
- swing-down dining table (5)
- side supports (6)
- restraining strap (6)
- padded foot support (6)
- wire basket (6)
- body supports (7)
- foot rest (7)



Merivaara Corp.

Puustellintie 2, FI-15150 LAHTI, FINLAND

Tel. +358 3 3394 611, fax +358 3 3394 6144

merivaara@merivaara.fi

www.merivaara.com