



The Gynecology/Urology Procedure Chair/Table

SONESTA 6303



The Stille-Sonesta 6303 chair... three motor positioning makes this a versatile, comfortable solution for any Gynecology and/or Urology procedure.



STILLE



Gynecology/Urology Procedure Chair
SONESTA

6303

The 6302 Fits Your Needs.

- General
- Examinations
- Cystoscopy
- Ultrasound
- Biofeedback
- Vasectomy
- OB/GYN Procedures
- Urodynamics



Engineered To Be Flexible.

- Three Motors for height, seat and back adjustments
- Easy positioning for patient and physician comfort
- Accepts the full range of Sonesta accessories
- Wide range of adjustment, from a low of 40cm (16") to 100cm (39.4") high
- Lightweight and mobile easy to move to other rooms; unmatched access for floor cleaning and disinfecting.
- Tough and reliable supports up to 375 lbs - 170 kg for years of service in the most demanding practices.

Standard Accessories:

- 3 motor control/foot plate 524-6357
- Cushion w/170mm cutout 524-6366-3
- Paper roll holder 524-6326
- 2 heavy duty Rissler leg rests 524-6336
- Drainage Receptacle 524-6325-1
- 2 hand grips/arm rests 524-6339-2
- 2 leg rest holders 524-6338

The Sonesta 6303 engineered to be the best, most flexible solution in your examination room.

The 6303 is an affordable path to Sonesta flexibility. Three motors provide full height, seat, and back adjustments. Packaged with basic accessories, the 6303 accepts the full line of Sonesta add-ons, allowing you to create the perfect table for your needs. The 6303's ergonomic design delivers a higher level of patient and physician comfort and flexibility... making it an incredible value in today's market!

SPECIFICATIONS:

LENGTH WITH EXTENSION CUSHION	175cm	69"
LENGTH WITHOUT	125cm	49.2"
WIDTH	65cm	25.5"
HEIGHT (MIN)	40cm	16"
HEIGHT (MAX)	100cm	39"
LIFT CAPACITY	170kg	375lbs
WHEEL DIAMETER (LOCKING)	10cm	4"
NET WEIGHT	<100kg	<220lbs

ELECTRICAL:

120/230v 24v DC Output
Meets EMC (Electromagnetic Compatibility) for IEC 60601
All Products carry the CE mark.



Stille Surgical AB
Sundbybergsvägen 1A
SE-171 73 Solna, Sweden
PHONE: +46 8 588 580 00
FAX: +46 8 588 580 05
Email: info@stille.se
www.stille.se

Distributed by: